

Columbus[®] AS Revision Knee System

with Advanced Surface Technology



Aesculap Orthopaedics

AESCULAP[®]
Implant Systems

Columbus[®] AS Revision Knee System

with Advanced Surface Technology



The patient's desire for a better quality of life has driven exponential demand in the knee replacement market and total knee revisions are projected to grow 601% between 2005 and 2030.¹ To address this increasing need, Aesculap Implant Systems offers the Columbus AS Revision Knee System that allows up to 130° of flexion and is based on the long standing technology of the Columbus Primary Knee System.


The Columbus AS Revision Knee System offers:

- Extensive portfolio of implant options for intraoperative flexibility
- Patented 7-layer Advanced Surface Technology:
 - provides substantial wear reduction, unmatched surface hardness and mechanical integrity.^{2,3}
 - designed for performance and potentially reduces metal ion release
- Precision engineering allows for high productivity and good wear rates
- Designed for use with the computer assisted navigation system, OrthoPilot[®] TKR



ADVANCED SURFACE TECHNOLOGY


Each layer of the proprietary 7-layer coating has a specific function. The outer layer, ZrN yields superior surface hardness and reduces abrasion. The 7-layer coating is specifically designed to bridge the difference in hardness from the ZrN over the CoCr substrate, improving the molecular structure of a layered system. This makes the molecular structure extremely stable against mechanical stresses and strains and results in a more resilient implant.^{11, 2}



The 5 transition layers provide stability and acts to create a barrier to help prevent metal ion release. The adhesive layer bonds the substrate securely to the coated layers, creating a firm connection between implant and subsequent coated layers.

Monolayer Coating

A hard surface on the relatively soft base material (CoCr) may lead to a higher risk of breakage of the surface, as it has been seen with monolayer coatings (eggshell effect).




Column structure of crystallines

VS

Multilayer Coating

The multilayer engineering in the transition layers leads to lower grain size and thus an improved elastic modulus, which specifically reduces shock to the material.



Small grain sizes

Columbus[®] AS Revision Knee System

with Advanced Surface Technology

PRECISION ENGINEERING

Variable offset can be adjusted to adapt to the patient's anatomy.

Advanced Surface Technology

- 7-layer coating delivers unparalleled surface hardness and greatly improved scratch resistance, as demonstrated in testing with other Aesculap knee designs.^{2,3}

Tibia offset

- Moves up to 6 mm offset in medial or lateral direction

Femur offset

- Moves up to 4 mm in anterior and posterior direction

Polyethylene Options


- MC (medium constrained) or varus/valgus stable HC (high constrained) are available.

4 Point locking mechanism

- Mitigates backside wear




BETA POLYETHYLENE DURABILITY



ASTM F 2003: artificial aging of 10 years acc. to Kurtz et al.⁶ 14 days/70°C/5bar O₂ shown in BiCONTACT study.

Improved age-resistance through beta radiation sterilization



	Gamma Sterilization	Beta Sterilization
Radiation	Lower intensity, deeper higher penetration, doses: 2.5 Mrad – 4 Mrad	Higher intensity, concentrated, lower penetration, doses: 2.5 Mrad – 4 Mrad
Sterilization Time	Longer: 16 hours	Shorter: 15 seconds
Result	Higher content of residual free radicals leading to a higher risk of oxidation	Fewer residual free radicals after sterilization process causing less oxidation ²⁷

Fig. 11: Gamma vs. Beta Sterilization

Columbus[®] AS Revision Knee System

with Advanced Surface Technology

INTRAOPERATIVE FLEXIBILITY

A wide variety of implants provide broad options for the patient's anatomy



2 polyethylene gliding surfaces

- Medium and High-Constrained
- Size range from 10 mm to 32 mm




4 patella sizes

- \varnothing 27 x 7
- \varnothing 30 x 8
- \varnothing 33 x 9
- \varnothing 36 x 10



7 femur sizes



femur spacers

- Interchangeable distal and posterior femur augments in 5, 10, 15 mm



6 tibia sizes

- Tibia hemi spacers in 5, 10, 15 mm
- Cemented and pressfit tibia stems in two lengths and seven sizes
 - Asymmetrical support flange for increased stability

\varnothing = diameter

FEMUR STEM OPTIONS


Asymmetrical support flange for increased stability




Cemented stem version



Pressfit stem version



4 longitudinal grooves to help avoid the risk of embolism



10 longitudinal grooves (Wagner Profile)

Femur extension stems

Cemented (6°): lengths 77 and 157 mm \varnothing 12, 15 and 18 mm
 Pressfit (5° / 7°): lengths 117 and 177 mm \varnothing 12, 13, 14, 15, 16, 17, 18, 19 and 20 mm

TIBIA STEM OPTIONS



Cemented stem version



Pressfit stem version

Tibia extension stems

Cemented: lengths 52 and 92 mm \varnothing 12, 15 and 18 mm
 Pressfit: lengths 92 and 132 mm \varnothing 11, 12, 13, 14, 15, 16, 17, 18, 19 and 20 mm

Indication:

The Columbus® Revision Knee System is indicated for use in reconstruction of the diseased knee joint caused by osteoarthritis, rheumatoid arthritis, post-traumatic arthritis, the need to revise failed arthroplasties or osteotomies where pain, deformity or dysfunction persist, and for patients suffering from correctable valgus or varus deformity and moderate flexion contracture. The Columbus Revision Knee System is designed for use with bone cement.

Contraindication:

Contraindications include, but are not limited to:

- Joint conditions that can be treated by reconstructive surgery (e.g. osteotomy)
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic diseases and metabolic disorders
- Severe osteoporosis or osteomalacia
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Bone malformations, axial misalignments or other bone conditions that rule out implantation of a prosthetic joint
- Predictable overload of the joint implant (e.g. due to adiposity)
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Fever, infection or inflammation (systemic or local)
- Pregnancy
- Mental illness
- Severe osteopenia (or any other medical or surgical finding) that would preclude any benefit from the implants
- Combination with implant components from other manufacturers
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials
- All cases listed under indications

References

1. J Bone Joint Surg Am. 2007 Apr 89(4): 780-5. Projections of primary and revision hip and knee arthroplasty in the United States from 2005 to 2030. Kurtz S, Ong K, Lau E, Mowat F, Halpern M.
2. Affatato S, Spinelli M, Lopomo N, Grupp TM, Marcacci M, Toni A. Can the method of fixation influence the wear behavior of ZrN coated unicompartmental mobile knee prostheses? Clin Biomech (Bristol, Avon). 2011 Feb; 26(2): 152-8. Epub 2010 Oct 8.
3. Grupp TM, Schwiesau T. Determination of the wear behavior of the UNIVATION mobile knee system T018, Mar 2007.
4. Bell CJ, Walker PS, Abeyundera MR, Simmons JM, King PM, Blunn GW. Effect of oxidation of delamination of ultrahighmolecular-weight polyethylene tibial components. J Arthroplasty 1998 Apr; 13(3):280-90.
5. Blömer W, Lohrmann E. Verschleißbeständigkeit von UHMWPE-Artikulationen in der Hüftgelenksendoprothetik. In: Weller S, Braun A, Eingartner C, Maurer F, Weise K, Winter E, Volkmann R. Das BICONTACT Hüftendoprothesensystem 1987-2007. Stuttgart: Georg Thieme Verlag; 2007. p. 94-100.
6. Sharkey PF, Hozack WJ, Rothman RH, Shastri S, Jacoby SM. Insall Award paper. Why are total knee arthroplasties failing today? Clin Orthop Relat Res. 2002 Nov;(404):7-13.
11. Aesculap data on file
25. Santana AE. Relating hardness-curve shapes with deformation mechanisms in TiAlN thin films indentation. Materials Science and Engineering A 406(2005) 11-18

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2014 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap Implant Systems, LLC | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 866-229-3002 | Fax 610-984-9096 | www.aesculapimplantsystems.com

Aesculap Implant Systems, LLC - a B. Braun company

DOC1047 Rev A 2M 9/14